

# DEVELOPMENT OF A RAPID HONEY SCREENING IMMUNOASSAY FOR RESIDUES OF NITROFURANS AOZ, AMOZ, AHD, SEM



Morris C., Fergie N., Kane N., Afrasiabi, Morteza.

Biorex Food Diagnostics, Unit 2C Antrim Technology Park, BT4 1 1QS, Antrim, Northern Ireland, United Kingdom

## BACKGROUND


Nitrofurans are a class of broad spectrum antibiotics which have potentially carcinogenic effects following human consumption. As a result, their use in food-producing animals has been banned in many jurisdictions including the European Union and USA. Reliable methods are therefore required to monitor their illegal use in food produce to ensure food safety.

Monitoring of the illegal use of Nitrofurans complicated by the short *in vivo* half life of the parent compounds. However, their tissue-bound metabolites – AMOZ, AOZ, AHD and SEM – are stable for several weeks after use of the parent compound and are therefore more reliable markers of the illegal use of Nitrofurans. This study reports the use of ELISA to quantify Nitrofuran metabolites in honey following their release from the sample matrix.

Parent Compound	Metabolite
Furazolidone	AOZ
Furaltadone	AMOZ
Nitrofurantoin	AHD
Nitrofurazone	SEM

## METHODS


Honey samples were prepared following a single sample preparation protocol for all four Nitrofuran metabolites, as illustrated in the diagram below.



Many of the reagents required for the sample preparation, including 2-NBA and spiking material, are provided in the assay kit. All reagents, including standards and antibodies required for the ELISA are also provided in ready-to-use form as part of the assay kit.

## RESULTS

**Calibration Range**  
**AOZ: 0.05 – 1.5 ppb**  
**AMOZ: 0.08 – 2.5 ppb**  
**AHD: 0.03 – 1.2 ppb**  
**SEM: 0.05 – 5.0 ppb**



**Figure 1:** Typical curves obtained using each Nitrofuran ELISA kit, with cross reactivities with other Nitrofuran metabolites highlighted in corresponding tables. The results highlight that cross-reactivities with other Nitrofuran metabolites are extremely low, enabling honey samples prepared in one run to be applied to all four Nitrofuran ELISA kits.

**AHD Intra-Assay Precision**

	Expected Concentration: 0.09ppb		Expected Concentration: 0.2ppb	
	ABSORBANCE (nm)	CONCENTRATION (ppb)	ABSORBANCE (nm)	CONCENTRATION (ppb)
MEAN	1.043	0.100	0.529	0.210
SD	0.04	0.01	0.02	0.01
%CV	4.2	6.7	4.3	3.5

**AMOZ Intra-Assay Precision**

	Expected Concentration: 0.08ppb		Expected Concentration: 0.25ppb	
	ABSORBANCE (nm)	CONCENTRATION (ppb)	ABSORBANCE (nm)	CONCENTRATION (ppb)
MEAN	1.778	0.100	0.938	0.340
SD	0.09	0.01	0.04	0.02
%CV	5.0	15.1	3.9	5.4

**AOZ Intra-Assay Precision**

	Expected Concentration: 0.125ppb		Expected Concentration: 0.5ppb	
	ABSORBANCE (nm)	CONCENTRATION (ppb)	ABSORBANCE (nm)	CONCENTRATION (ppb)
MEAN	1.883	0.122	0.455	0.714
SD	0.09	0.01	0.03	0.04
%CV	4.6	9.0	6.7	5.6

**SEM Intra-Assay Precision**

	Expected Concentration of 0.05ppb		Expected Concentration: 0.5ppb	
	ABSORBANCE (nm)	CONCENTRATION (ppb)	ABSORBANCE (nm)	CONCENTRATION (ppb)
MEAN	1.887	0.081	0.422	0.800
SD	0.06	0.01	0.04	0.12
%CV	3.4	10.7	9.5	14.6

**Figure 2:** Intra-assay precision for each of the four Nitrofurans. The results demonstrate good intra-assay precision for all four Nitrofuran ELISAs, with CV for absorbance all much less than 10%

## RESULTS

AHD Inter-Assay Precision			AMOZ Inter-Assay Precision		
Expected Concentration: 0.2 ppb			Expected Concentration: 0.1 ppb		
MEAN	0.200		MEAN	0.090	
SD	0.02		SD	0.01	
%CV	9.6		%CV	10.0	

AOZ Inter-Assay Precision			SEM Inter-Assay Precision		
Expected Concentration: 0.25 ppb			Expected Concentration: 0.25 ppb		
MEAN	0.250		MEAN	0.270	
SD	0.02		SD	0.00	
%CV	9.9		%CV	9.3	

**Figure 3:** Inter-assay precision for each of the four Nitrofurans. Inter-assay precision was examined by running one control level across up to twenty independent test runs. Inter-assay precision is good across all four Nitrofuran kits, with %CV ≤10%.

## CONCLUSIONS

Biorex Food Diagnostics have developed four ELISA kits for the detection of four individual Nitrofuran Metabolites. Each of the kits exhibits very low cross-reactivity with the other Nitrofuran metabolites, enabling one sample preparation to be used across all four kits.

Each of the Nitrofuran ELISA kits have a number of advantages over other commercially-available kits:

- **Short assay time:** 45 minutes
- **Shorter sample derivatisation time** of 90 minutes as opposed to 120 minutes
- **One sample preparation** can be used across all four Nitrofuran kits
- **Reduced solvent volumes.** This offers a number of benefits including:
  - Lower reagent costs
  - Faster sample evaporation
  - Preparation can be completed in one 15ml tube
- Spiking material is provided in the assay kit in a **ready-to-use** format
- Conjugate is supplied **ready-to-use**, removing the need for dilution prior to use
- The **derivatisation reagent** 2-nitrobenzaldehyde is supplied as a **ready-to-use** liquid